

250100 LEA Symbols Folding Pediatric Eye Chart - Answer Key

Testing distance of 10ft/3m

Distance Equivalents
 10 Foot 3 Meter
 Letter Size

20 Foot 6Meter
 Decimal LogMar

10/100 30M	3/30											20/200 .10	6/60 1.0
10/80 24M	3/24											20/160 .12	6/48 0.9
10/63 19M	3/19											20/125 .16	6/38 0.8
10/50 15M	3/15											20/100 .20	6/30 0.7
10/40 12M	3/12											20/80 .25	6/24 0.6
10/32 9.5M	3/9.5											20/63 .32	6/19 0.5
10/25 7.5M	3/7.5											20/50 .40	6/15 0.4
10/20 6M	3/6											20/40 .50	6/12 0.3
10/16 4.8M	3/4.8											20/32 .63	6/9.5 0.2
10/12.5 3.8M	3/3.8											20/25 .80	6/7.5 0.1
10/10 3M	3/3											20/20 1.0	6/6 0.0
10/8 2.4M	3/2.4											20/16 1.25	6/4.8 -0.1
10/6.3 1.9M	3/1.9											20/12.5 1.6	6/3.8 -0.2
10/5 1.5M	3/1.5											20/10 2.0	6/3 -0.3
10/4 1.2M	3/1.2											20/8 2.5	6/2.4 -0.4