

500072 - Landolt C Distance Vision Answer Key

Testing distance of 10ft/3m

Actual Size								Equivalent			
10 Foot	3 Meter	Letter Size						20 Foot	6 Meter	Decimal	LogMAR
10/100	3/30	30M	⦶	⦶	⦶	⦶	⦶	20/200	6/60	.10	1.0
10/80	3/24	24M	⦶	⦶	⦶	⦶	⦶	20/160	6/48	.12	0.9
10/63	3/19	19M	⦶	⦶	⦶	⦶	⦶	20/125	6/38	.16	0.8
10/50	3/15	15M	⦶	⦶	⦶	⦶	⦶	20/100	6/30	.20	0.7
10/40	3/12	12M	⦶	⦶	⦶	⦶	⦶	20/80	6/24	.25	0.6
10/32	3/9.5	9.5M	⦶	⦶	⦶	⦶	⦶	20/63	6/19	.32	0.5
10/25	3/7.5	7.5M	⦶	⦶	⦶	⦶	⦶	20/50	6/15	.40	0.4
10/20	3/6	6M	⦶	⦶	⦶	⦶	⦶	20/40	6/12	.50	0.3
10/16	3/4.8	4.8M	⦶	⦶	⦶	⦶	⦶	20/32	6/9.5	.63	0.2
10/12.5	3/3.8	3.8M	⦶	⦶	⦶	⦶	⦶	20/25	6/7.5	.80	0.1
10/10	3/3	3M	⦶	⦶	⦶	⦶	⦶	20/20	6/6	1.0	0.0
10/8	3/2.4	2.4M	⦶	⦶	⦶	⦶	⦶	20/16	6/4.8	1.25	-0.1
10/6.3	3/1.9	1.9M	⦶	⦶	⦶	⦶	⦶	20/12.5	6/3.8	1.6	-0.2
10/5	3/1.5	1.5M	⦶	⦶	⦶	⦶	⦶	20/10	6/3	2.0	-0.3