

# 756300 - Near Vision Answer Key

## Testing distance of 16in/40cm

### Distance Equivalents

Foot	Meter		Letter Size	Decimal
20/400	6/120	R N O V S	8.0M	.05
20/320	6/95	Z C R D H	6.3M	.063
20/250	6/75	N V S O K	5.0M	.08
20/200	6/60	D R Z K O	4.0M	.10
20/160	6/48	S N H C V	3.2M	.12
20/125	6/38	C R V S Z	2.5M	.16
20/100	6/30	V K C N H	2.0M	.20
20/80	6/24	S V K D N	1.6M	.25
20/63	6/19	K D H Z C	1.25M	.32
20/50	6/15	H Z C O R	1.0M	.40
20/40	6/12	O K D H N	.80M	.50
20/32	6/9.5	Z O N K C	.63M	.63
20/25	6/7.5	R H S V D	.50M	.80
20/20	6/6	D S O R Z	.40M	1.0
20/16	6/4.8	H Z S C V	.32M	1.25
20/12.5	6/3.8	C K R D Z	.25M	1.6
20/10	6/3	R D O N K	.20M	2.0

# 756300 - Near Vision Answer Key

## Testing distance of 16in/40cm

### Distance Equivalents

Foot	Meter		Letter Size	Decimal
20/400	6/120	D S R K N	8.0M	.05
20/320	6/95	C K Z O H	6.3M	.063
20/250	6/75	O N R K D	5.0M	.08
20/200	6/60	K Z V D C	4.0M	.10
20/160	6/48	V S H Z O	3.2M	.12
20/125	6/38	H D K C R	2.5M	.16
20/100	6/30	C S R H N	2.0M	.20
20/80	6/24	S V Z D K	1.6M	.25
20/63	6/19	N C V O Z	1.25M	.32
20/50	6/15	R H S D V	1.0M	.40
20/40	6/12	S N R O H	.80M	.50
20/32	6/9.5	O D H K R	.63M	.63
20/25	6/7.5	Z K C S N	.50M	.80
20/20	6/6	C R H D V	.40M	1.0
20/16	6/4.8	N V K D O	.32M	1.25
20/12.5	6/3.8	K Z D C R	.25M	1.6
20/10	6/3	D N R K O	.20M	2.0