

# 904034 - Distance Vision Answer Key

## Testing distance of 10ft/6m

### Distance Equivalentents

Foot	Meter		Letter Size	Decimal
20/40	6/12	○ □ ○ ♥ ♠	.80M	.50
20/40	6/12	♠ ○ ♥ □ ♥	.80M	.50
20/32	6/9.5	□ ○ ♠ ♥ □	.63M	.63
20/32	6/9.5	○ ♥ ♠ □ ○	.63M	.63
20/25	6/7.5	♠ ○ ♥ ♠ □	.50M	.80
20/25	6/7.5	○ □ ♥ ♠ ○	.50M	.80
20/20	6/6	♠ ♥ □ ○ ♥	.40M	1.0
20/20	6/6	□ ♥ ♠ □ ○	.40M	1.0

# 904034 - Near Vision Answer Key

## Testing distance of 16in/40cm

### Distance Equivalents

Foot	Meter							Letter Size	Decimal					
20/200	6/60					□	△	♥	○	△	4.0M	.10		
20/40	6/12	○	□	○	△	♥		♥	△	□	○	△	.80M	.50
20/32	6/9.5	□	○	□	♥	△		□	♥	△	□	○	.63M	.63
20/50	6/15	♥	○	♥	△	□		♥	○	□	♥	△	1.0M	.40
20/40	6/12	□	♥	○	□	△		○	□	○	△	♥	.80M	.50
20/32	6/9.5	○	△	□	♥	○		□	○	□	♥	△	.63M	.63
20/25	6/7.5	△	♥	△	○	□		○	□	♥	△	○	.50M	.80
20/20	6/6	□	△	○	♥	△		△	♥	□	○	♥	.40M	1.0
20/16	6/4.8	○	□	♥	△	♥		♥	△	○	□	○	.32M	1.25