



Brock String

Technical Bulletin

Traditional Brock Strings with plastic Beads



Brock Strings consist of a white nylon string with five or three small plastic beads of different colors. They are used for treatment of convergence insufficiency and other anomalies of binocular vision as well as to disrupt suppression. During therapy, the end with the stick is held near the tip of the nose while the other end with the loop is hung to a fixed point like a door knob. Dozen pack includes 12 sets of instructions.

Length	Number of Beads	Individual P/N	Twelve Pack P/N
20 Ft	3	603100	593700
12 Ft	5	603300	602400
10 Ft	3	603200	593800
6 Ft	5	603400	593900



Brock String

Technical Bulletin

Cord Lock Brock Strings - stay in place



Beads that fail to stay in place are the bane of traditional Brock Strings. The cord stretches or the exercise requires a sloped string and then the beads start moving. Our new Cord Lock Brock Strings solve the problem by using a spring-loaded bead that is guaranteed to stay in place! During therapy, the end with the stick is held near the tip of the nose while the other end with the loop is hung to a fixed point like a door knob. Our new handle looks sharp and helps the user focus on the task, not the string. Just like the traditional Brock Strings, Cord Lock versions are used for treatment of convergence insufficiency and other anomalies of binocular vision as well as to disrupt suppression. All of the standard configurations are offered. Each pack comes with instructions for each Brock String.

Length	Number of Beads	Individual P/N	3 Pack P/N	Twelve Pack P/N
20 Ft	3	702801		702812
12 Ft	5	703001		703012
10 Ft	Depends on P/N	702901 (with 5 beads)	603500 (with 5 beads)	702912 (with 3 beads)
6 Ft	5	703101	703103	703112