

# Projection Mirrors, Screens and Mounts



## Swivel Mount Mirrors

### Front Surface Mirror Kit

Permits the use of smaller rooms to accomplish a full testing distance with a projector. The front surface maximizes reflected illumination and minimizes distortion. 705xx Kit includes two 12" (31cm) square mirrors. Each mirror includes a universal ball joint mount. Kit is based on projection onto painted wall and does not include a screen. 6702xx is a single mirror from the 705x x kit.

**705xx Front Surface Kit**

**6702xx Single Front Surface Mirror**



Swivel Mount Mirrors support projection of 20 ft (6M) charts in rooms shorter than 20 feet (6M). These Rear coated Mirrors are 12 inch (30 cm) square. Black metal frame and Swivel Mount are sturdy for many years of use and blend with any décor. Swivel mount permits extensive flexibility in mounting at different angles in both axis directions.

**5693R Single Swivel Mirror Rear coated**

**5746R Pair of Swivel Mirror Rear coated**



Unframed Screen

## Projection Screens And Mounts

Non-polarizing unframed aluminum screens. The wall or swivel mounts can be fastened to the back of the screens using a two-part epoxy glue (not included). Screens are washable. The assembly can be mounted to the ceiling or a wall.

**4452R Unframed Screen 13" (34 cm) Sq.**

**4453R Unframed Screen 16" (41 cm) Sq.**

**4454R Unframed Screen 20" (51 cm) Sq.**



Swivel Mount

**4455R Swivel Mount**

**4456R Black Wall Mount**



Black Wall Mount